

High Performance Fume Hood

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Aerodynamic Sash Pull	48"	Air	Resin-Chem*	Yes	3"x6"	Up to 16"	Yes	Yes
Streamlined Airflow Design	60"	Gas						
Combination Horizontal/ Vertical Sash Pull	72"	Cold Water						
Locking Sash Stop	84"	Vacuum						
Secondary Containment Trough	96"							
Air Foil and Air Flow Monitor								
2 ea. Std. 115v Duplex Outlet								

**Optional liner materials are available depending on chemical uses.*

Bypass RAV (reduced air volume) Fume Hood

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Reduces Sash Opening to Reduce Air Volume	48"	Air	Resin*	Yes	3"x6"	Up to 16"	Yes	Yes
Face Velocity 60-150	60"	Gas						
Interior Lighting/Light Switch	72"	Cold Water						
Air Flow Monitor	84"	Vacuum						
2 ea. Std. 115v Duplex Outlet	96"							

**Optional liner materials are available depending on chemical uses.*

Perchloric Acid Fume Hood*

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Wash Down Spray Bar	48"	Air	Type 316 Coved Stainless Steel	Yes	No	No	Yes**	Yes
Face Velocity 60-150	60"	Gas						
Integral Work Top	72"	Cold Water						
Integral Trough in Rear	84"	Vacuum						
Interior Lighting/Light Switch	96"							
2 ea. Std. 115v Duplex Outlet								

**Specifically designed for use with Perchloric Acid*

***Standard and digital alarm types available*

Radioisotope Acid Fume Hood*

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Welded Seams	48"	Air	Type 316 Coved Stainless Steel	Yes	No	No	Yes**	Yes
Face Velocity 60-150	60"	Gas						
Integral Work Top	72"	Cold Water						
Water Tight Surface that Supports 200lbs./sq.in.	84"	Vacuum						
Interior Lighting/Light Switch	96"							
2 ea. Std. 115v Duplex Outlet								

**Specifically designed for use with radioisotopes*

***Standard and digital alarm types available*

Walk-In Fume Hood

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Able to Roll in Large Amounts of Equipment	48"	Air	Resin-Chem	Yes	3"x6"	No	Yes**	Yes
Interior Height of 633/4"	60"	Gas						
Accommodate Complex Set-ups	72"	Cold Water						
Two Independently Raising Vertical Sashes	84"	Vacuum						
Interior Lighting/ Light Switch	96"							
2 ea. Std. 115v Duplex Outlet								

***Standard and digital alarm types available*

Add-Air Fume Hood

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Built-in Bypass to Maintain Velocity	48"	Air	Resin-Chem	Yes	3"x6"	No	No	No
Provides up to 70% Overall Air Volume	60"	Gas						
Face Velocity 60-150	72"	Cold Water						
	84"	Vacuum						
	96"							

Variable Air Volume (VAV) Fume Hood

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Front Sash Opening	48"	Air	Resin-Chem	Yes	3"x6"	No	No	No
Lintel Restricted Air Intake	60"	Gas						
Face Velocity 60-150	72"	Cold Water						
One Finger Sash Operation	84"	Vacuum						
Exhaust Volume Varied by VAV Controls	96"							

Demonstration Fume Hood

Features	Sizes	Remote Fixtures	Hood Liners	Fan	Cup Sink	Sink	Alarm	Custom Sizes
Sash Incorporated on Both Sides	48"		Resin-Chem	Yes	3"x6"	No	Yes**	No
End Viewing Window	60"							
Demonstration Can be Viewed from 3 Sides	72"							
	84"							
	96"							

***Standard and digital alarm types available*

Cabinets

Type	Sizes (inches)
Acid Cabinets	18", 24", 30", 36", 42", 48"
Solvent Cabinets	18", 24", 30", 36", 42", 48"
Standard Cabinets	18", 24", 30", 36", 42", 48"